


# Colourimeter

# Tainlab



- Three LED sources – Red, Blue and Green
- No need to change filters or replace globes
- Gives direct reading of absorbance
- Determine concentrations using Beer's Law
- Monitor progress of reactions
- Print graphs and save results to file
- Reload calibrations from results files
- Includes PC software and 20 cuvettes
- Works with any Tainlab interface

**Tain**  
ELECTRONICS  
Established 1982

10 Rowern Court, Box Hill North, Vic 3129, Australia  
Tel: 03 9898 7366 Fax: 03 9899 7992  
email: [info@tain.com.au](mailto:info@tain.com.au)  
web: [www.tainlab.com](http://www.tainlab.com)  
ABN: 70 006 496 706